## InvivoDental 7.2

US: August 13<sup>th</sup>, 2025

<https://learn.osteoidinc.com/invivo-user-manuals>

Anatamage, Inc.  
3350 Scott Blvd. Bldg. #29  
Santa Clara, CA 95054, U.S.A.  
Tel: +1 (408) 333-3484  
Email: support@osteoidinc.com  
www.osteoidinc.com

Freyr Regulatory GmbH  
Marie-Curie-Straße 8,  
D-79539 Lörrach, Germany  
Tel: +49 618 170 79007  
www.freyrsolutions.com

**EU:**  
Freyr Regulatory GmbH  
Marie-Curie-Straße 8,  
D-79539 Lörrach, Germany  
Tel: +49 618 170 79007  
EU-Importer@freyrsolutions.com  
www.freyrsolutions.com

**Australia Sponsor  
(TGA)**

Henry Schein Halas  
Locked Bag 5003  
Alexandria, NSW, 2015  
customer.care@henryschein.com.au  
Tel: +61 1300-65-88-22

**South Korea**

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<p>Indications for Use</p>	<p>InvivoDental is a software application used for the display and 3D visualization of medical image files from scanning devices, such as CT, MRI, or 3D Ultrasound. It is intended for use by radiologists, clinicians, referring physicians, and other qualified individuals to retrieve, process, render, review, store, print, assist in diagnosis and distribute images, utilizing standard PC hardware. Additionally, InvivoDental is a preoperative software application used for the simulation and evaluation of dental implants, orthodontic planning, and surgical treatments.</p> <p>This device is not indicated for mammography use.</p>
<p>Minimum System Requirements</p>	<p>Windows 11</p> <ul style="list-style-type: none"> <li>• CPU Processor: 64-bit, Intel i3 or AMD Ryzen 3</li> <li>• RAM Memory: 8 GB</li> <li>• GPU Graphics: Dedicated Nvidia/AMD Graphics, 1GB Memory</li> <li>• Storage: 100GB</li> <li>• Display Resolution Scale: Not to exceed 150%</li> </ul> <p>MacOS 13 (Ventura, 64-bit)</p> <ul style="list-style-type: none"> <li>• 2019 Hardware and Later</li> </ul>
<div style="text-align: center;">  </div> <p>Safety Warnings and Precautions – Refer to ANA003 Invivo 7.2 Instructions for Use</p>	<ul style="list-style-type: none"> <li>• The software provides tools for dental prosthetics creation and design but is dependent on the user to determine and use the appropriate parameters. Incorrect parameters may affect the quality of the final milled products or otherwise to surgical delays or complications.</li> <li>• Please ensure all hardware devices with the software are password-protected from unauthorized use and all patient information is secure.</li> <li>• Saving files in lossy format degrades image quality.</li> <li>• InVivoDental’s measurement functionality has been verified to be accurate up to +/- 0.10 mm. However, measurement accuracy depends on the image data and the hardware scanner that generated the image data. The measurement cannot be more precise than the resolution of the image. Software reports the value based on user-picked points. Due to the nature of medical imaging, the boundary is not always well defined. The apparent boundary depends on the current brightness and contrast setting. The boundary may shift as the user makes adjustments to brightness and contrast. The user must understand the limitation of the measurement value before applying to the patient. Any measurement that is incorrect can lead to surgical complications if diagnosis, treatment plans and/or actual treatment is based on the incorrect measurements. It is critical for the end user to learn how to perform measurements correctly and employ proper usage of all measurement tools. If you notice any inconsistencies or software problems with measurements, or have further questions or concerns about using measurement tools correctly, please contact us at (408) 333-3484 or email us at support@osteoidinc.com.</li> <li>• The measurement values may not be true anatomical volumetric measurements. Due to the nature of the imaging, there are imaging artifacts such as white noise, scattering, beam hardening, ring noise or off scale H.U. The software measurement tool cannot distinguish the imaging artifact from the true anatomy. Furthermore, the measurement values depend on threshold values; thus, user must set proper thresholds to get the best estimate of the desired structure. The measurements must not be used as the sole metric for any treatment.</li> <li>• Any nerve that is traced in a way that does not conform to the actual pathway of the nerve can lead to surgical complications if diagnosis, treatment plans and/or actual treatment is based off of the incorrect tracing. It is critical for the end user to learn how to properly perform nerve tracings correctly. If you notice any inconsistencies or software problems with nerve tracing or have further questions or concerns about nerve tracing, please contact Osteoid support at (408) 333-3484 or email us at support@osteoidinc.com</li> </ul>

- Any implants and restorations that are incorrectly visualized in Invivo may lead to surgical complications or treatment delays. If you notice any inconsistencies or software problems with implant or restoration planning or have further questions or concerns, please contact Osteoid support at (408) 333-3484 or email us at support@osteoidinc.com
- Any implant planning that is performed incorrectly in any aspect including but not limited to implant location, orientation, angle, diameter, length, and/or manufacture can lead to surgical complications if diagnosis, treatment plans and/or actual treatment is based off of the implant planning error(s) in question. It is critical for the end user to learn how to use the implant treatment planning tools correctly. If you notice any inconsistencies or software problems with implant planning or have further questions or concerns about correct utilization of implant planning, please contact Osteoid support at (408) 333-3484 or email us at support@osteoidinc.com
- Incorrect crown fitting may cause inaccurate crown shape and positioning. This may lead to surgical complications or delayed if treatment planning.
- Sequence direction depends on the focal trough direction. Medial and lateral may be flipped if the focal trough angle is in the opposite direction.
- Soft tissue prediction is an approximation and may not reflect actual results.
- 3DAnalysis module is not a medical device and should be used as a reference or presentation tool only.
- Results produced by the software's diagnostic and planning tools are dependent on correct selection of landmarks by the user in the 3D volume. The software reports values based on these user-picked points. Due to the nature of medical imaging, the boundaries are not always well defined. The apparent boundaries depend on the current brightness and contrast settings. The boundaries may shift as the user makes adjustments to brightness and contrast. The user must understand the limitations of the measurement values before applying to the patient. If you notice any inconsistencies or software problems with measurements, or have further questions or concerns about correctly using measurement tools correctly, please contact us at (408) 333-3484 or email us at support@osteoidinc.com
- The analysis results are dependent on the norm values that are configured within each configuration file, which may contain predefined as well as user-defined data. The predefined norms that are available upon installation may or may not be based on accepted clinical values for the respective measurements. Users must verify the validity of norm values that will be used before applying this function on patient data. For further questions regarding the function, configuration, and validity of norm data, please contact the Osteoid Support Department at (408) 333-3484 or support@osteoidinc.com.
- 3D Surgery tool's intended use is for patient consultation and statistical analysis only. For diagnoses, use Invivo.
- If the control is created outside of the page borders, an error message is displayed, and the control is not created.
- When resizing a Slice Group Control, note that the slice numbers may disappear if the control is downsized past a limit (approximately the size that the numbers can no longer fit within each respective slice).
- Bitmaps under 24 bit may not display correctly.
- Stitching module is not a medical device and should be used as a reference or presentation tool only
- Treatment information can only be manually saved. All information added after last saving will be lost.

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